



## Stock Glass Diagnostic Screw Neck Vials

Tubular glass vials manufactured by SCHOTT from SCHOTT Fiolax® neutral Type I glass tubing. Vials can be manufactured to your precise specification for viable production runs




Tubular Type I glass diagnostic/screw neck vials

### CLEAR VIALS

Adelphi Part No.	Size	a	b	c	e	m	Pack Size
VCD003	3ml	18	36	14	8.6	22.3	273
VCD005	5ml	18	41.5	14	8.6	22.3	285
F-3059	7ml	22	40	18	11.5	24	186
VCD010	10ml	23.75	48	18	11.5	37	160
VCD020S18	20ml	30	55	18	11.5		91

### AMBER VIALS

Adelphi Part No.	Size	a	b	c	e	m	Pack Size
VAD003	3ml	18	36	14	8.6	<b>22.3</b>	285
VAD005	5ml	18	41.5	14	8.6	<b>22.3</b>	285
1009754	7ml	22	40	18	11.6	<b>23.2</b>	186
VADDIN010	10ml	24	45	18	11.5	<b>29</b>	154
VAD010S22	10ml	28	55	22	15.2	<b>35</b>	112
VA24X60S18AT	15ml	24	60	18	11.5		154
VAD030S22	30ml	33.5	67	22	15.2	<b>46</b>	80



All measurements in (mm) and volumes nominal

### SCHOTT Quality

- All production steps are monitored by a DIN EN ISO 9001 GMP oriented Quality Management System
- All Manufacturing processes are carried out in accordance with ISO 15378
- All vials are subjected to opto-electronic dimension inspection during the production process
- Packing is carried out in class 100,000 clean room conditions
- In addition, SCHOTT vials are supported by an FDA Drug Master File (DMF) Type III registration

Vials are available siliconised for viable quantities